

Tru Products



Tru Abutment (CAD/CAM)

TruAbutment CAD/CAM custom abutments provide individual anatomic contours, emergence profile and margin placement to provide an excellent foundation for improved long-term function and esthetics compared to stock abutments.

- **Patient-specific**

CAD/CAM abutments are individually designed, taking the final tooth shape into consideration.


- **Adjustable cement line**

Better maintenance of proper biological width.

- **High-quality, cost-effective**

Our manufacturing process makes sure you get a precise fit every time.

** Screws for CAD/CAM and Titanium Base are not inter-changeable*



Tru Base

Designed Ideally for Zirconia Restoration
Available in engaging and non-engaging

** for CAD/CAM Use Only*

Tru Scan Body

Hybrid: PEEK+Titanium
Can be used intraorally
Usable in autoclave



Tru Lab Analog

Stainless Steel
Non-anodized



Compatible with

Astra Tech™ TX



System (Color)	3.0 ●	3.5/4.0 ●	4.5/5.0 ●
Library Code	Tru_Astra 3.0	Tru_Astra 3.5/4.0	Tru_Astra 4.5/5.0
Screw Torque	15N.Cm	20N.Cm	25N.Cm
Screw Size	M1.4 x 0.3P	M1.6 x 0.35P	M2.0 x 0.4P
Screw Driver	Hex 1.27mm (0.05 inch)		
Tru Abutment (CAD/CAM)	AS3-CS	AS4-CS	AS5-CS
Tru Ti-Base	Single - Short	AS3-B13S	AS4-B13S
	Bridge - Short	AS3-B13B	AS4-B13B
	Single - Long	AS3-B15S	AS4-B15S
	Bridge - Long	AS3-B15B	AS4-B15B
Tru Scan Body	AS3-SB	AS4-SB	AS5-SB
Tru Lab Analog	AS3-LA	AS4-LA	AS5-LA

Compatible with

Biomet 3i™ Certain



System (Color)	3.4 ●	4.1 ●	5.0 ●	6.0 ●
Library Code	Tru_3i Certain 3.4	Tru_3i Certain 4.1	Tru_3i Certain 5.0	Tru_3i Certain 6.0
Screw Torque	20N.Cm			
Screw Size	M1.6 x 0.35P			
Screw Driver	Hex 1.2mm (0.048 inch)			
Tru Abutment (CAD/CAM)	BC34-CS	BC41-CS	BC50-CS	BC60-CS
Tru Ti-Base	Single - Short	BC34-B13S	BC41-B13S	BC50-B13S
	Single - Long	BC34-B15S	BC41-B15S	BC50-B15S
Tru Scan Body	BC34-SB	BC41-SB	BC50-SB	BC60-SB
Tru Lab Analog	BC34-LA	BC41-LA	BC50-LA	BC60-LA

Compatible with

Dentis

OneQ-SL (S-Clean)




System (Color)	Narrow ●	Regular ●
Library Code	Tru_Dentis OneQ-SL (Sub) Mini	Tru_Dentis OneQ-SL (S-Clean) Regular
Screw Torque	15N.Cm	30N.Cm
Screw Size	M1.4 x 0.3P	M2.0 x 0.4P
Screw Driver	Hex 1.27mm (0.05 inch)	
Tru Abutment (CAD/CAM)	DQN-CS	DQR-CS
Tru Ti-Base	Single - Short	DQN-B13S
	Bridge - Short	DQN-B13B
	Single - Long	DQN-B15S
	Bridge - Long	DQN-B15B
Tru Scan Body	DQN-SB	DQR-SB
Tru Lab Analog	DQN-LA	DQR-LA

Compatible with

Dentis

OneQ-SL (I-Clean)



System (Color)	RN ●	WN ●
Library Code	Tru_Dentis OneQ-SL (i-Clean RN)	Tru_Dentis OneQ-SL (i-Clean WN)
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	DIR-CS	DIW-CS
Tru Ti-Base	Single - Short	DIR-B13S
	Bridge - Short	DIR-B13B
	Single - Long	DIR-B15S
	Bridge - Long	DIR-B15B
Tru Scan Body	DIR-SB	DIW-SB
Tru Lab Analog	DIR-LA	DIW-LA

Compatible with


Dentium SuperLine



System (Color)	SuperLine ●	
Library Code	Tru_Dentium SuperLine	
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.27mm (0.05 inch)	
Tru Abutment (CAD/CAM)	DSU-CS	
Tru Ti-Base	Single - Short	DSU-B13S
	Bridge - Short	DSU-B13B
	Single - Long	DSU-B15S
	Bridge - Long	DSU-B15B
Tru Scan Body	DSU-SB	
Tru Lab Analog	DSU-LA	

Compatible with

Dentium SimpleLine II



System (Color)	RN ●	WN ●
Library Code	Tru_Dentium SimpleLine II RN	Tru_Dentium SimpleLine II WN
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.27mm (0.05 inch)	
Tru Abutment (CAD/CAM)	DSR-CS	DSW-CS
Tru Ti-Base	Single - Short	DSR-B13S
	Bridge - Short	DSR-B13B
	Single - Long	DSR-B15S
	Bridge - Long	DSR-B15B
Tru Scan Body	DSR-SB	DSW-SB
Tru Lab Analog	DSR-LA	DSW-LA

Compatible with

DIO UF



System (Color)	Mini (Narrow Fix) ●	Regular ●
Library Code	Tru_Dio UF Mini	Tru_Dio UF Regular
Screw Torque	20N.Cm	30N.Cm
Screw Size	M1.6 x 0.35P	M2.0 x 0.4P
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	DUN-CS	DUR-CS
Tru Ti-Base	Single - Short	DUN-B13S
	Bridge - Short	DUN-B13B
	Single - Long	DUN-B15S
	Bridge - Long	DUN-B15B
Tru Scan Body	DUN-SB	DUR-SB
Tru Lab Analog	DUN-LA	DUR-LA

Compatible with

DIO Internal



System (Color)	RN ●	WN ●
Library Code	Tru_Dio internal RN	Tru_Dio Internal WN
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	DRN-CS	DWN-CS
Tru Ti-Base	Single - Short	DRN-B13S
	Bridge - Short	DRN-B13B
	Single - Long	DRN-B15S
	Bridge - Long	DRN-B15B
Tru Scan Body	DRN-SB	DWN-SB
Tru Lab Analog	DRN-LA	DWN-LA

Compatible with

DIO SM



System (Color)	Narrow ●	Regular ●
Library Code	Tru_Dio SM Narrow	Tru_Dio SM Regular
Screw Torque	35N.Cm	
Screw Size	M1.8 x 0.35P	M2.0 x 0.4P
Screw Driver	Torx 1.7mm	
Tru Abutment (CAD/CAM)	DMS-CS	DMR-CS
Tru Ti-Base	Single - Short	DMS-B13S
	Bridge - Short	DMS-B13B
	Single - Long	DMS-B15S
	Bridge - Long	DMS-B15B
Tru Scan Body	DMS-SB	DMR-SB
Tru Lab Analog	DMS-LA	DMR-LA

Compatible with

Friadent (Xive®)



System (Color)	3.0 ●	3.4 ●	3.8 ●	4.5 ●	5.5 ●
Library Code	Tru_Friadent(Xive) 3.0	Tru_Friadent(Xive) 3.4	Tru_Friadent(Xive) 3.8	Tru_Friadent(Xive) 4.5	Tru_Friadent(Xive) 5.5
Screw Torque	24N.Cm				
Screw Size	M1.4 x 0.3P	M1.6 x 0.35P			
Screw Driver	Hex 1.2mm (0.048 inch)				
Tru Abutment (CAD/CAM)	DX30-CS	DX34-CS	DX38-CS	DX45-CS	DX55-CS
Tru Ti-Base	Single - Short	DX30-B13S	DX34-B13S	DX38-B13S	DX45-B13S
	Single - Long	DX30-B15S	DX34-B15S	DX38-B15S	DX45-B15S
Tru Scan Body	DX30-SB	DX34-SB	DX38-SB	DX45-SB	DX55-SB
Tru Lab Analog	DX30-LA	DX34-LA	DX38-LA	DX45-LA	DX55-LA

Compatible with

Hiossen HG/ET

(Osstem GS/TS)




System (Color)	Mini ●	Regular ●
Library Code	Tru_Hiossen ET Mini (Osstem TS Mini)	Tru_Hiossen ET Standard (Osstem TS Standard)
Screw Torque	20N.Cm	30N.Cm
Screw Size	M1.6 x 0.35P	M2.0 x 0.4P
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	OTM-CS	OTR-CS
Tru Ti-Base	Single - Short	OTM-B13S
	Bridge - Short	OTM-B13B
	Single - Long	OTM-B15S
	Bridge - Long	OTM-B15B
Tru Scan Body	OTM-SB	OTR-SB
Tru Lab Analog	OTM-LA	OTR-LA

Compatible with

Hiossen SS

(Osstem SS)



System (Color)	RN ●	WN ●
Library Code	Tru_Hiossen SS RN (Osstem SS RN)	Tru_Hiossen SS WN (Osstem SS WN)
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	OSR-CS	OSW-CS
Tru Ti-Base	Single - Short	OSR-B13S
	Bridge - Short	OSR-B13B
	Single - Long	OSR-B15S
	Bridge - Long	OSR-B15B
Tru Scan Body	OSR-SB	OSW-SB
Tru Lab Analog	OSR-LA	OSW-LA

Compatible with

Keystone Dental Prima



System (Color)	SD ●	RD ●	WD ●
Library Code	Tru_Keystone Prima SD	Tru_Keystone Prima RD	Tru_Keystone Prima WD
Screw Torque	30N.Cm		
Screw Size	M1.8 x 0.35P		
Screw Driver	Square Driver		
Tru Abutment (CAD/CAM)	KPS-CS	KPR-CS	KPW-CS
Tru Ti-Base	Single - Short	KPR-B13S	KPW-B13S
	Single - Long	KPR-B15S	KPW-B15S
Tru Scan Body	KPS-SB	KPR-SB	KPW-SB
Tru Lab Analog	KPS-LA	KPR-LA	KPW-LA

Compatible with

Megagen Anyone



System (Color)	Anyone ●	
Library Code	Tru_Megagen Anyone	
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	MAR-CS	
Tru Ti-Base	Single - Short	MAR-B13S
	Bridge - Short	MAR-B13B
	Single - Long	MAR-B15S
	Bridge - Long	MAR-B15B
Tru Scan Body	MAR-SB	
Tru Lab Analog	MAR-LA	

Compatible with

Megagen Anyridge



System (Color)		Anyridge ●
Library Code		Tru_Megagen Anyridge
Screw Torque		30N.Cm
Screw Size		M1.8 x 0.35P
Screw Driver		Hex 1.2mm (0.048 inch)
Tru Abutment (CAD/CAM)		MRD-CS
Tru Ti-Base	Single - Short	MRD-B13S
	Bridge - Short	MRD-B13B
	Single - Long	MRD-B15S
	Bridge - Long	MRD-B15B
Tru Scan Body		MRD-SB
Tru Lab Analog		MRD-LA

Compatible with

Megagen Exfeel



System (Color)		Exfeel ●
Library Code		Tru_Megagen Exfeel
Screw Torque		30N.Cm
Screw Size		M2.0 x 0.4P
Screw Driver		Hex 1.2mm (0.048 inch)
Tru Abutment (CAD/CAM)		MEF-CS
Tru Ti-Base	Single - Short	MEF-B13S
	Bridge - Short	MEF-B13B
	Single - Long	MEF-B15S
	Bridge - Long	MEF-B15B
Tru Scan Body		MEF-SB
Tru Lab Analog		MEF-LA

Compatible with

MIS Seven (M4)



System (Color)	Narrow ●	Standard ●	Wide ●
Library Code	Tru_MIS Seven (M4) Narrow	Tru_MIS Seven (M4) Standard	Tru_MIS Seven (M4) Wide
Screw Torque	30N.Cm		
Screw Size	M1.6 x 0.35P	1-72UNF	
Screw Driver	Hex 1.27mm (0.05 inch)		
Tru Abutment (CAD/CAM)	MIN-CS	MIS-CS	MIW-CS
Tru Ti-Base	Single - Short	MIS-B13S	MIW-B13S
	Single - Long	MIS-B15S	MIW-B15S
Tru Scan Body	MIN-SB	MIS-SB	MIW-SB
Tru Lab Analog	MIN-LA	MIS-LA	MIW-LA

Compatible with


Neobiotech IS



System (Color)	IS ●	
Library Code	Tru_Neobiotech IS	
Screw Torque	30N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Hex 1.2mm (0.048 inch)	
Tru Abutment (CAD/CAM)	NIS-CS	
Tru Ti-Base	Single - Short	NIS-B13S
	Bridge - Short	NIS-B13B
	Single - Long	NIS-B15S
	Bridge - Long	NIS-B15B
Tru Scan Body	NIS-SB	
Tru Lab Analog	NIS-LA	

Compatible with

Neobiotech IT



System (Color)	RN ●	WN ●	
Library Code	Tru_Neobiotech IT RN	Tru_Neobiotech IT WN	
Screw Torque	30N.Cm		
Screw Size	M2.0 x 0.4P		
Screw Driver	Hex 1.2mm (0.048 inch)		
Tru Abutment (CAD/CAM)	NTR-CS	NTW-CS	
Tru Ti-Base	Single - Short	NTR-B13S	NTW-B13S
	Bridge - Short	NTR-B13B	NTW-B13B
	Single - Long	NTR-B15S	NTW-B15S
	Bridge - Long	NTR-B15B	NTW-B15B
Tru Scan Body	NTR-SB	NTW-SB	
Tru Lab Analog	NTR-LA	NTW-LA	

Compatible with

Nobel Biocare Active®



System (Color)	3.0 ●	NP ●	RP ●	WP ●	
Library Code	Tru_Nobel Active 3.0	Tru_Nobel Active NP	Tru_Nobel Active RP	Tru_Nobel Active WP	
Screw Torque	15N.Cm	35N.Cm			
Screw Size	M1.4 x 0.3P	M1.6 x 0.35P	M2.0 x 0.4P		
Screw Driver	Unigrip				
Tru Abutment (CAD/CAM)	NA3-CS	NAN-CS	NAR-CS	NAW-CS	
Tru Ti-Base	Single - Short	NA3-B13S	NAN-B13S	NAR-B13S	NAW-B13S
	Bridge - Short	NA3-B13B	NAN-B13B	NAR-B13B	NAW-B13B
	Single - Long	NA3-B15S	NAN-B15S	NAR-B15S	NAW-B15S
	Bridge - Long	NA3-B15B	NAN-B15B	NAR-B15B	NAW-B15B
Tru Scan Body	NA3-SB	NAN-SB	NAR-SB	NAW-SB	
Tru Lab Analog	NA3-LA	NAN-LA	NAR-LA	NAW-LA	

Compatible with

Nobel Biocare Replace®



Tru Abutment
(CAD/CAM)



Tru Ti-Base



Tru Scan Body



Tru Lab Analog

System (Color)	NP ●	RP ●	WP ●	60 ●
Library Code	Tru_Nobel Replace NP	Tru_Nobel Replace RP	Tru_Nobel Replace WP	Tru_Nobel Replace 60
Screw Torque	35N.Cm			
Screw Size	M1.8 x 0.35P	M2.0 x 0.4P		
Screw Driver	Unigrip			
Tru Abutment (CAD/CAM)	NRN-CS	NRR-CS	NRW-CS	NR6-CS
Tru Ti-Base	Single - Short	NRN-B13S	NRW-B13S	NR6-B13S
	Single - Long	NRN-B15S	NRW-B15S	NR6-B15S
Tru Scan Body	NRN-SB	NRR-SB	NRW-SB	NR6-SB
Tru Lab Analog	NRN-LA	NRR-LA	NRW-LA	NR6-LA

Compatible with

Straumann® Bone Level



Tru Abutment
(CAD/CAM)



Tru Ti-Base



Tru Scan Body



Tru Lab Analog

System (Color)	NC ●	RC ●
Library Code	Tru_Straumann Bone Level NC	Tru_Straumann Bone Level RC
Screw Torque	35N.Cm	
Screw Size	M1.6 x 0.35P	
Screw Driver	Torx	
Tru Abutment (CAD/CAM)	SBN-CS	SBR-CS
Tru Ti-Base	Single - Short	SBN-B13S
	Bridge - Short	SBN-B13B
	Single - Long	SBN-B15S
	Bridge - Long	SBN-B15B
Tru Scan Body	SBN-SB	SBR-SB
Tru Lab Analog	SBN-LA	SBR-LA

Compatible with

Straumann® Tissue Level



System (Color)	RN ●	WN ●
Library Code	Tru_Straumann Tissue Level RN	Tru_Straumann Tissue Level WN
Screw Torque	35N.Cm	
Screw Size	M2.0 x 0.4P	
Screw Driver	Torx	
Tru Abutment (CAD/CAM)	STR-CS	STW-CS
Tru Ti-Base	Single - Short	STR-B13S
	Bridge - Short	STR-B13B
	Single - Long	STR-B15S
	Bridge - Long	STR-B15B
Tru Scan Body	STR-SB	STW-SB
Tru Lab Analog	STR-LA	STW-LA

Compatible with

Zimmer® Screw Vent



System (Color)	3.1 ●	3.5 ●	4.5 ●	5.7 ●
Library Code	Tru_Zimmer Screw Vent 3.1	Tru_Zimmer Screw Vent 3.5	Tru_Zimmer Screw Vent 4.5	Tru_Zimmer Screw Vent 5.7
Screw Torque	30N.Cm			
Screw Size	M1.6 x 0.35P	1-72UNF		
Screw Driver	Hex 1.27mm (0.05 inch)			
Tru Abutment (CAD/CAM)	ZVN-CS	ZV3-CS	ZV4-CS	ZV5-CS
Tru Ti-Base	Single - Short	ZVN-B13S	ZV3-B13S	ZV4-B13S
	Bridge - Short	ZVN-B13B	ZV3-B13B	ZV4-B13B
	Single - Long	ZVN-B15S	ZV3-B15S	ZV4-B15S
	Bridge - Long	ZVN-B15B	ZV3-B15B	ZV4-B15B
Tru Scan Body	ZVN-SB	ZV3-SB	ZV4-SB	ZV5-SB
Tru Lab Analog	ZVN-LA	ZV3-LA	ZV4-LA	ZV5-LA